

# Pulseq on GE

Jon-Fredrik Nielsen

November 16, 2023

# Outline

- 1 Overview
- 2 Representing the .seq file using 'segments'
- 3 Write segment representation to file
- 4 Simulate sequence in Pulse Studio
- 5 Execute on scanner
- 6 Misc topics: 'runs' feature; TTL trigger

---

## Pulseq on GE v1.0 (DRAFT)

Instructions for running Pulseq sequences on GE  
using the TOPPE universal pulse sequence interpreter

This document applies to version 6 of the interpreter (tv6.e).

*This version of the manual contains GE specific information  
and can be distributed to sites with a GE research scanner.*

This document was last modified on Sun Nov 5 16:32:56 2023

LaTeX source for this manual can be found at:  
<https://github.com/jfnrielsen/TOPPEpdSourceCode/tree/UserGuide/v6/>

An occasionally updated PDF version of this manual can be found at:  
<https://drive.google.com/file/d/1eToTYtFip6lUaAoh0furoldF14h8cDs/view?usp=sharing>

---

- Redacted version for the public: <https://toppemri.github.io/>
- Full version in the MR Software Sharing forum on [gecares.com](https://gecares.com)

# History

Original 'TOPPE' interpreter <sup>1</sup>:

- Developed independently of Pulseseq; tailored to GE
- Proved that cross-vendor sequence programming is possible
- Shortcomings:
  - ▶ gap between blocks ( $\simeq 200\text{-}400\mu\text{s}$ )
  - ▶ Gradients had to start and end at 0 within each block
  - ▶ Different version of interpreter source code for each scanner software version
- Workarounds possible

---

<sup>1</sup>Nielsen JF, Noll DC. TOPPE: A framework for rapid prototyping of MR pulse sequences. Magnetic resonance in medicine. 2018 Jun;79(6):3128-34.

# Pulseq philosophy: assemble blocks on the fly


- Define blocks on the fly during sequence execution (streaming)

```
# Format of blocks:  
# NUM DUR RF  GX  GY  GZ  ADC  EXT  
[BLOCKS]  
1 316    1   0   0   1   0   1  
2 100    0   2   3   4   0   0  
3 223    0   0   0   0   0   0  
4  46    0   5   0   0   1   0  
5 100    0   6   7   8   0   0  
6 256    0   0   0   0   0   0
```

- Not directly compatible with EPIC


# New Pulseseq interpreter (2023)

- Time gaps eliminated (mostly)
- Gradient waveforms can cross block boundaries
- A single set of source files for multiple software versions
- Reuses some basic TOPPE functionality (e.g., file formats), but mostly rewritten from scratch
- Sequence represented as a collection of scaled 'segments'



# Segment concept

- Segment = a fixed sub-sequence of Pulseseq blocks that are always executed together
- Can consist of any number of blocks (even 100s)



# Convert a .seq file to segment representation

- Step 1: Label the start of each segment<sup>2</sup>

```
segmentID = 1; % any non-negative integer
seq.addBlock(rf, gz, mr.makeLabel('SET', 'LIN', segmentID));
```

- Step 2: Convert

```
>> ceq = seq2ceq('gre2d.seq');
ceq =
      nMax: 192
nParentBlocks: 4
parentBlocks: {[1x1 struct] [1x1 struct] [1x1 struct] [1x1 struct]}
      segments: [1x1 struct]
      nSegments: 1
      loop: [192x10 double]
```

'ceq' object:


- A vendor-independent sequence description
- May find use in other contexts (e.g., other vendor platforms, simulation tools, etc)

---

<sup>2</sup>Need to do manually for now



```
>> seq.plot('blockRange', [1 8], 'showBlocks', true);
```



# 'Parent' blocks

Unique Pulseseq blocks (up to a scaling/phase factor)

```
>> ceq.parentBlocks{1}
ans =

    struct with fields:

        blockDuration: 0.0032
            rf: [1x1 struct]
            gx: []
            gy: []
            gz: [1x1 struct]
            adc: []
        label: [1x1 struct]
        amp: [1x1 struct]
```

# Segment definitions

In our example, there's only one segment:

```
>> ceq.segments
```

```
ans =  
  
    struct with fields:  
        groupID: 1  
        nBlocksInGroup: 6  
        blockIDs: [1 2 0 3 4 0]    % parent block IDs (0 = wait period)
```

# Scan loop definition

Dynamic scan information such as waveform amplitude scaling, RF/ADC phase, etc

```
>> ceq.loop(1:8,:)

```

```
ans =
```

```
1.0e+06 *
```

0.0000	0.0000	0.0000	0	0	0	0	0.2667	0
0.0000								
0.0000	0.0000	0	0	0	-0.1105	-0.0907	-0.4312	0
0.0000								
0.0000	0	0	0	0	0	0	0	0
0.0000								
0.0000	0.0000	0	0	0	0.5556	0	0	0
0.0000								
0.0000	0.0000	0	0	0	1.6931	0.0907	1.7021	0
0.0000								
0.0000	0	0	0	0	0	0	0	0
0.0000								
0.0000	0.0000	0.0000	0.0000	0	0	0	0.2667	0
0.0000								
0.0000	0.0000	0	0	0	-0.1105	-0.0850	-0.4312	0
0.0000								

# Write segment representation to file

- Convert to GE hardware units and write to file:

```
>> ceq2ge(ceq, sysGE, 'gre2d.tar');
```

```
ceq2ge: Writing block 192/192  
preflightcheck: Checking max 10s SAR, time interval 0–10s (scan duration: 0s)  
    Predicted peak 10s SAR in 150 lbs subject: 0.0 W/kg  
Sequence file gre2d.tar ready for execution on GE scanners
```

- Both steps can be combined into one:

```
>> seq2ge('gre2d.seq', sysGE, 'gre2d.tar');
```

- The .tar file can now be simulated and executed on GE scanners

# Summary of file creation steps

## 1 Define the start of each segment in the .seq file

```
segmentID = 1;    % any non-negative integer  
seq.addBlock(rf, gz, mr.makeLabel('SET', 'LIN', segmentID));
```

## 2 Convert to .tar file

```
>> seq2ge('grd2d.seq', sysGE, 'gre2d.tar');
```

# Simulate sequence in Pulse Studio

## Pulse Studio

- GE's new sequence simulator, available for software version MR30.2
- Use WTools to create waveform files, then view with Pulse View
- Display whole (or parts of) sequence: waveforms, timing, SSP pulses
- GE distributes the Python source for Pulse View
  - ▶ Wants to collaborate with the community to improve Pulse View

## Simulation steps

- 1 Compile interpreter for simulation (MR30.2)
- 2 Simulate scan loop in WTools (instructions in Pulseseq on GE manual)
- 3 Load waveform files in Pulse View

# Simulate sequence in Pulse Studio





## Execute on scanner

- 1 Copy .tar file to folder of your choice (/any/path/)
- 2 Place toppe<ID>.entry file in [/usr/g/research/pulseq/v6/](#) (<ID> = integer of your choice)
- 3 Replace the 1st line in toppe<ID>.entry with /any/path/
- 4 Set sequence control variables (CVs): TODO
- 5 Scan

Steps 1-4 are the same as for the Pulse Studio simulator

Additional instructions in Pulseq on GE manual

# Pulseq on GE workflow summary

## 1 Define the start of each segment in the .seq file

```
segmentID = 1; % any non-negative integer  
seq.addBlock(rf, gz, mr.makeLabel('SET', 'LIN', segmentID));
```

## 2 Convert to .tar file

```
>> seq2ge('grd2d.seq', sysGE, 'gre2d.tar');
```

## 3 Simulate in Pulse Studio

- 1 Compile interpreter for simulation (MR30.2)
- 2 Simulate scan loop in WTools
- 3 Load waveform files in Pulse View

## 4 Run sequence on scanner

- 1 Compile interpreter for your scanner software version
- 2 Prescribe and run

# Dynamic imaging using 'runs' feature

- Dynamic sequences (e.g., fMRI) can be very long → large .seq files
- Redundancy: typically a short sub-sequence is simply repeated
  - ▶ But pay attention to RF phase cycling (e.g., RF spoiling)
- Pulseseq exploits this on the interpreter level
  - ▶ The .seq file only needs to define the sub-sequence
  - ▶ Set number of repetitions ('runs') on UI
  - ▶ Shorter scan file → loads faster on scanner
  - ▶ New feature currently in 'develop' branch

# TTL trigger output pulses

- TTL output pulses are specified in the Pulseseq file as follows:

```
trig_out = mr.makeDigitalOutputPulse('ext1', 'delay', 100e-6);  
seq.addBlock(gx, trig_out);    % block containing a gradient and a TTL trigger
```

- seq2ge.m captures TTL triggers
- New GE interpreter:
  - ▶ Trigger<sup>3</sup> starts at the requested time
  - ▶ Trigger duration defaults to 500us, use 'trigdur' control variable (CV) to change
  - ▶ Trigger channel (DABOUT) is set on UI
  - ▶ Trigger cannot be the only event inside a block (may change in future)
  - ▶ Built-in support for Skope calibration scan<sup>4</sup>

---

<sup>3</sup>Implemented by Dinank Gupta

<sup>4</sup>Thanks to Andy Derbyshire for this feature

# Acknowledgments

## Interpreter code and testing:

- Rolf Schulte (GE) provided C code for PNS calculation, advised on implementation of safety checks, and provided example code for dynamic waveform loading (via the fidall psd which he is the author of).
- Kang Wang (GE) consulted on various implementation and usability aspects, e.g., FOV shifting and C++ code integration.
- Ante Zhu and Afis Ajala (GE) helped with testing and debugging on the MAGNUS head gradient system.
- Sherry Huang (GE) has been central to connecting the development team with core GE resources.
- Dinank Gupta added TTL trigger out support.
- John Derbyshire contributed code for adding the Skope field camera calibration scan in MPS2.
- Matteo Cencini added son-of-recon CV and example scripts.
- Dirk Poot added error handling code to an early version (v2).
- Alexei Ouriadov added broadband support in an early version (v2)

## MATLAB toolboxes:

- Amos Cao wrote the convenient `toppe.write2loop()` MATLAB function, and made many other contributions.
- Melissa Haskell made a number of contributions, including model-based reconstruction support for spiral imaging.
- Rolf Schulte (GE) provided the `toppe.pns()` MATLAB script for PNS calculation.
- Github user moedn proposed use of the MIT license

If you feel someone should be added to this list, please let [jfnielse@umich.edu](mailto:jfnielse@umich.edu) know – any omissions are entirely unintentional!