

# ICE and Gadgetron Reconstruction Based on the Pulseseq Framework

**Qingping Chen**

*Division of Medical Physics, Dept. Of Radiology,  
University Medical Center Freiburg, Germany*

*Nov 16, 2023*

# Outline

- Workflow for data acquisition and image reconstruction
- How to enable ICE/Gadgetron (based on `s21_GRE2D_GRAPPA`:  
[https://github.com/pulseseq/ISMRM-Virtual-Meeting--November-15-17-2023/tree/main/tutorials/day2\\_imageReconstruction](https://github.com/pulseseq/ISMRM-Virtual-Meeting--November-15-17-2023/tree/main/tutorials/day2_imageReconstruction) )
  - ADC labeling & label debugging
  - sequence definitions
  - interface setting
- Example reconstructions
  - MPRAGE with GRAPPA = 2
  - Multi-slice EPI with navigator and ramp sampling

# Workflow


- Sequence design in the Pulseseq software:
  - ADC labeling
  - sequence definitions (`seq.setDefinition(...)`)
- Extended Pulseseq interpreter *sequence*
  - Load ADC labels and sequence definitions
  - Link to ICE and Gadgetron<sup>1,2</sup>
- ICE is integrated to Siemens, while Gadgetron is vendor-independent by employing ISMRMRD<sup>3</sup> data format (Siemens, GE, Philip, Bruker, etc).

Gadgetron: <https://github.com/gadgetron/gadgetron>

ISMRMRD: [https://github.com/ismrmrd/ge\\_to\\_ismrmrd](https://github.com/ismrmrd/ge_to_ismrmrd)

IceGadgetron: [https://github.com/NHLBI-MR/ISMRM2019\\_demo/tree/master/Siemens\\_Scanner\\_Setup](https://github.com/NHLBI-MR/ISMRM2019_demo/tree/master/Siemens_Scanner_Setup)


1. Hansen, et al., MRM, 2013. 2. Xue, et al., MRM, 2015. 3. Inati, et al. MRM, 2017.



# ICE/Gadgetron Setup

## ADC labeling


- Phase-encoded **index** labeling (caution: c-style numbering!):
  - Index setting
    - `mr.makeLabel( 'SET', 'LIN', 0 )`
  - Index increment
    - `mr.makeLabel('INC', 'LIN', 2 )`
- ACS **flag** labeling:
  - Reference & image scan
    - `mr.makeLabel('SET', 'REF', true )`
    - `mr.makeLabel('SET', 'IMA', true )`
  - Reference scan
    - `mr.makeLabel('SET', 'REF', true )`
    - `mr.makeLabel('SET', 'IMA', false )`
  - Image scan
    - `mr.makeLabel('SET', 'REF', false )`
    - `mr.makeLabel('SET', 'IMA', false )`



# ICE/Gadgetron Setup

## ADC debugging

```
adc_lbl =  
seq.evalLabels('evolution','adc');  
figure ;  
hold on ;  
plot(adc_lbl.LIN) ;  
plot(adc_lbl.REF) ;  
plot(adc_lbl.IMA) ;
```




LIN: PE index

# ICE/Gadgetron Setup

## ADC debugging

```
adc_lbl =  
seq.evalLabels('evolution','adc');  
figure ;  
hold on ;  
plot(adc_lbl.LIN) ;  
plot(adc_lbl.REF) ;  
plot(adc_lbl.IMA) ;
```



REF: ACS reference scan flag

IMA: ACS reference&image scan flag



# ICE/Gadgetron Setup

## Sequence definitions

- **2D multi-slice mode**

```
seq.setDefinition('SlicePositions', 0) ;  
seq.setDefinition('SliceThickness', sliceThickness) ;  
seq.setDefinition('SliceGap', sliceGap) ;
```


- **Additional information required by GRAPPA**

```
seq.setDefinition('kSpaceCenterLine', centerLineIdx-1) ;  
seq.setDefinition('PhaseResolution', phaseResoluion) ;
```

# ICE/Gadgetron Setup

## Interface setting

- Turn on iPAT setting: *PAT mode*
- Set *Accel. factor PE*, *Ref. lines PE*, and *Reference scan mode*
- **Caution:** GRAPPA mode must be turned on for ICE/Gadgetron recon. **Do not** forget to turn it off after data acquisition.






# ICE/Gadgetron Setup

## Interface setting


- Select “Data handling”
  - none: no data handling
  - auto export: NIH data catcher
  - ICE 2D: 2D recon (VB, VE, etc)
  - ICE 3D: 3D recon (VB, VE, etc)
  - ICE STD: 2D&3D recon (XA, etc)
  - MRD-export: ISMRMRD<sup>1</sup> data export
  - gadget2D: gadgetron online recon
- Caution: IceGadgetron<sup>2</sup> needs to be installed in the Siemens scanner for online Gadgetron reconstruction.



1. Inati, et al. MRM, 2017.; 2. Xue, et al., MRM, 2015.


# Example reconstructions

MPRAGE with GRAPPA = 2




# Example reconstructions

MPRAGE with GRAPPA = 2



# Example reconstructions


## Multi-slice EPI with navigator and ramp sampling





# Example reconstructions

## Multi-slice EPI with navigator and ramp sampling



## Acknowledgements:

Berkin Bilgic

Frank Zijlstra

Jon-Fredrik Nielsen

Moritz Zaiss

Qiang Liu

Sebastian Littin

Borjan Gagoski

Imam Shaik

Juergen Hennig

Naveen Murthy

Maxim Zaitsev

Will Grissom

Douglas Noll

Jeff Fessler


Mojtaba Shafiekhani

Niklas Wehkamp

Scott Peltier

Yogesh Rathi

# THANK YOU FOR YOUR ATTENTION!



Supported by NIH U24-NS120056 (Nielsen, Zaitsev)  
and R01-EB032378 (Rathi, Bilgic)