

What's New Pulseq

Maxim Zaitsev

*Division of Medical Physics, Dept. of Radiology,
University Medical Center Freiburg, Germany*

Nov 16, 2023

What is 'New' in a Continuum?

- Some interesting/relevant developments from the past ~12 months
- RF simulation
- Peripheral nerve stimulation (PNS)
- Data labels, slice positions
- Mechanical resonances
- Some more fun stuff...

RF Simulation & Analysis

In addition to TOPPE's slice profile calculation...

mr.calcRfCenter() & **mr.calcRfBandwidth()**



- return RF peak position / FT width
- Fast, only for approximate guidance, small tip-angle approximation (STA)

[Mz_z, Mz_xy, F, ref_eff, Mx_xy, My_xy]=

mr.simRf(rfPulse, preph_fact, reph_fac)

- Fast RF simulation based on the rotation formalism (ignoring relaxation)

See [demoUnsorted/demoRfSimulation.m](#) in the main Pulseq repository for numerous use examples




Peripheral Nerve Stimulation Prediction

- Direct integration with the SAFE model (Siemens-specific)
- Open-SAFE by Filip Szczepankiewicz
https://github.com/filip-szczepankiewicz/safe_pns_prediction

seq.calcPNS()

- Analyze the sequence in memory for PNS
- Ideally requires access to the pertinent scanner-specific PNS model data (stored in *.asc files)



Pulseq Extensions

- Pulseq file format readily supports very flexible Extension objects
- Established extensions
 - Cardiac triggering
 - Trigger (digital output) pulses
 - Data labeling (the traditional use of the LABEL extension)
- Labels to control sequence execution
 - ONCE labels
 - NOPOS, NOROT, NOSCALE labels
- Labels may provide hints for the interpreters
 - New GE (and possibly candidate Philips) interpreters detect the “TR Loop”

New Label Display and Debugging Options



- Labeling of the raw data is often desirable
 - Online recon & integration with Gadgetron
- `seq.plot('TimeDisp', 'ms', 'label', 'lin,slc');`
- `adc_lbl=seq.evalLabels('evolution','adc');`

`adc_lbl =`

struct with fields:

```

LIN: [0 0 2 4 6 8 10 12 14 16 18 20
IMA: [0 0 0 0 0 0 0 0 0 0 0 1 0 1 0
REF: [0 0 0 0 0 0 0 0 0 0 0 1 1 1 1
NOISE: [1 0 0 0 0 0 0 0 0 0 0 0 0 0 0
SLC: [0 0 0 0 0 0 0 0 0 0 0 0 0 0 0
    
```




```

7 38 39 40 41 42 43 44 46 48 50 52 54 56 58 60 ... ]
0 0 0 0 0 0 0 0 0 0 0 0 0 1 0 1 0 1 0 1 0 ... ]
0 0 0 0 0 0 0 0 0 0 0 0 0 1 1 1 1 1 1 1 1 ... ]
0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 ... ]
0 0 0 0 0 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 ... ]
    
```

New Labels: ONCE


- Some interpreters offer the option run/repeat the Pulseseq sequence multiple times without a delay
- Label ONCE when set to
 1. Mark blocks that are only executed on the first repetition
 2. Mark blocks that are only executed on the last repetition
- Possible use cases:
 - Calibration scans performed once at the beginning
 - Signal preparation
 - dummy scans in GRE or
 - $\alpha/2$ TR/2 preparation in trueFISP
- As example see [demoSeq/writeTrufi.m](#)




Multi-Slice Mode

- Traditionally Pulseseq presents itself as a “3D sequence”
 - FOV in the [Definitions] section inscribes all slices
 - FOV positioning works similarly for all sequences
 - Positions of individual slices are not visible to the user
- New definitions switch the sequence to a 2D multi-slice mode:
seq.setDefinition('SlicePositions', slicePositionsArray);
seq.setDefinition('SliceThickness', thickness_in_m);
seq.setDefinition('SliceGap', sliceGap_in_m);

3D-mode vs. 2D-mode



'legacy' 3D mode




new 2D multi-slice mode

- Homogeneous handling for all imaging and spectroscopy sequences
- Inner life of the sequence is opaque

- 2D sequences need special handling
- Correct slice position information in raw data
- More familiar user experience

Gradient Spectrogram Analysis

- CAUTION! Gradient coil have **mechanical resonances** (forbidden echo spacings)
- Pulseseq now provides a demo (not a function yet) for a spectral analysis of the sequence in memory, see [demoUnsorted/gradSpectrum.m](#)
- We would like to hear from the manufacturers how to handle this knowledge





Ready for the Fun Part?

Pulseq & Gradient Music

- Experimental **mrMusic** library for creating gradient music

DISCLAIMERS:

- **You** are **completely, entirely and absolutely** responsible for the health or damage of your gradient system
- **Always** check the gradient spectrum and **avoid** the **mechanical resonances** of your system
- **Never** use the **maximum** gradient **amplitude**

mrMusic how-to

- Located in demoUnsorted of the main repository (lib + 4 demo scripts)
 - musicBadinerie.m, musicEntertainer.m, musicRootBeerRag.m & musicSoundTest.m
- 1. Initialize the mrMusic environment:
`mrMusic.init; % this creates over 30 global variables`
- 2. Create the melody as a cell array of bars and channels (n x 3)
`melody = { [...], [...], [...] ; }`
- 3. Convert melody to a channel-frequency table
`[pitches, durations] = mrMusic.melodyToPitchesAndDurations(...);`
- 4. Populate the Sequence object with the music piece
`seq = mrMusic.musicToSequence(seq, ...);`
- 5. Check how it sounds
`seq.sound();`

More about Melody and Notes

```
melody = {  
    % bar 1  
    [ o/4 h2/8 d3/16 h2/16 ],...  
    [ o/4 d2/8 o/8 ],...  
    [ o/4 hbb/8 hbb/8 ];...  
    % bar 2  
    [ fis2/8 h2/16 fis2/16 d2/8 fis2/16 d2/16 ], ...  
    [ h1/8 o/8 fis1/8 o/8 ], ...  
    [ db/8 hbb/8 fisb/8 db/8 ]; }
```

- Notes: a, b, c, d, e, f, g, h; ‘_is’ is high ‘_es’ is flat; ‘a1’ is the middle A
- Octaves: _bb, _b, _1, _2, _3;
- ‘o’ is a delay
- To define duration divide the note by the fractional duration, e.g. cis2/16 is C# of the second octave with the duration of 1/16th
- **Enjoy!**

mrMusic: Concluding Remarks

- Really easy to create nice 3-voice tunes
- Music can be combined with imaging or spectroscopy
- Many more options available (explore examples)
 - Sound quality settings (wave vs. sawtooth)
 - Temperaments for Baroque music
 - Staccato & Legato
 - **Tuning for avoiding mechanical resonances**
 - ...
- **Enjoy responsibly!**

Acknowledgements:

Berkin Bilgic

Frank Zijlstra

Jon-Fredrik Nielsen

Moritz Zaiss

Qiang Liu

Sebastian Littin

Borjan Gagoski

Imam Shaik

Juergen Hennig

Naveen Murthy

Qingping Chen

Will Grissom

Douglas Noll

Jeff Fessler

Mojtaba Shafiekhani

Niklas Wehkamp

Scott Peltier

Yogesh Rathi

THANK YOU FOR YOUR ATTENTION!



Supported by NIH U24-NS120056 (Nielsen, Zaitsev)
and R01-EB032378 (Rathi, Bilgic)

